

# VICTOR®

## Model T9000

### Non-Contact Infrared Forehead & Wrist Thermometer

Victor's T9000 durable non-contact infrared thermometer is quick, hygienic and convenient to use. It is perfect for all settings, including home, work/office, school, and hospital environments. This contactless thermometer reads forehead and wrist temperatures from 2-3 inches away, avoiding any direct physical contact with skin and possible cross-infection between multiple people. The T9000 can be used to read temperatures of adults and kids. The touchless and precise sensor makes it the perfect temperature reader, even on infants and toddlers while they are asleep.



#### Color Alert Display

FDA Approved



#### FEATURES

- **Contactless Technology:** Reads temperature from forehead or wrist with no physical contact (2-3 in.) avoiding direct contact with skin and possible cross-infection between multiple people.
- **Fast Readings:** 1-second readings.
- **Accurate Readings:** Measurement accuracy is within 0.4 degrees Fahrenheit.
- **Dual Modes (Human and Object):** In addition to measuring Human forehead and wrist temperature, the T9000 can switch to Object Mode and measure object, room, and liquid temperatures with the high precision infrared sensor.
- **32 Memory Storage:** View 32 previous readings; press "Mem" button continuously to display last 32 readings.
- **F/C Modes:** Switch between Fahrenheit and Celsius.
- **3-Color Backlit Display:**
  - When temperature is between 93.2 – 99.3 degrees Fahrenheit, backlit is Green
  - When temperature is between 99.5 – 100.2 degrees Fahrenheit, backlit is Yellow
  - When temperature is between 100.4 – 109.4 degrees Fahrenheit, backlit is Red
- **Sound and Color Alert:** Notifies when high temperature/fever detected.
- **Easy Trigger Operation:** Point and press. Beeper will sound for Completed Scan.
- **Battery Level Indicator on Display:** When battery is low, indicator light will flash. (Comes with 2 AAA)
- **Durable:** Thermometer has met drop test and continues to work properly even after accidental drop.
- **Automatic Power Off:** After 25 seconds of non-use.

*\*Default factory settings turn on to Human and Fahrenheit modes.*

[www.victortech.com](http://www.victortech.com) (800) 628-2420